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
The screenshot displays the spectrometer control interface. At the top, the 'Jobリスト' (Job List) shows a tree view with 'sample (1)' expanded to show 'Proton' with a duration of 0:02. A red circle highlights the 'Experiment' icon in the top toolbar. An arrow points from this icon to the 'Experiment Parameters' dialog box, which is also highlighted with a red circle. The dialog box contains various settings for the experiment, including storage filename, filename, comment, and acquisition parameters. A red 'X' icon in the top right corner of the dialog indicates the close button. Below the dialog, the '測定Queue' (Measurement Queue) is visible, showing a table with columns for 'サンプル' (Sample) and 'Status'. The status bar at the bottom provides system information such as laser gain, rotation speed, lock status, temperature, and liquid levels.

Header	Instrument	Acquisition	Pulse	ダイアグラム	Favorites
storage_filename	save_fname	save_fname			
filename	proton				
comment	single_pulse				
auto_filter	<input checked="" type="checkbox"/>				
auto_gain	<input type="checkbox"/>				
filter_limit	8				
force_tune	<input type="checkbox"/>				
save_aborted	<input checked="" type="checkbox"/>				
fsqd_factor	1				

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Info
Gain Value Established: -8
Starting Experiment:
Starting Collection
Building Output File: 256 kB
Sending file to data server
Post-experiment Default Initialization

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